Produktbeschreibung

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

Achtung

Product description

FotoCast® is a radiation-curing resin based on methacrylate for the generative manufacturing of forms with the Digital Light Processing (DLP®) technology for producing earmolds. For processing FotoCast® a DLP® system and a high performance light curing unit like e.g. PCU LED or Sonolux Flash Evo with nitrogen flushing is needed.

Processing

- To prepare a construction job, the parameters for FotoCast® have to be adjusted in the specific software appropriate for the DLP® system that is used.
- Check the filling of the material in the DLP® system, before starting a construction. If necessary, fill up the resin reservoir of the DLP® system with FotoCast® to the manufacturer's recommended level. We recommend wearing rubber gloves during the application of FotoCast®.
- After the construction job is completed and the platform has reached its parking position, a dripping off time of approx. 1 minute is recommended.
- We recommend cleaning the FotoCast® shell with Isopropanol for approx. 4 minutes. After that the FotoCast® shell should be dried with compressed air.
- Post-curing
 - PCU LED 30 sec. / 50% light output, incl. vacuum
 - Sonolux Flash Evo 5 minutes gas pre-flow + 20 sec. curing
- First the inside of the FotoCast® shell should be injected with separator and secondly dried with compressed air.
- Like this, the manufactured FotoCast® shell can optimally be used for the forming of e.g. Biopor® silicone products.
- After the silicone is completely cured, it can be easily removed from the FotoCast® shell by peeling the shell off manually, without any additional tools.
- Due to wrong application for instance, a fracture of the component in the construction process might occur that could also lead to a contamination of the material. However, FotoCast® can still be filtered. Furthermore, it is recommended to clean the resin reservoir of the DLP® system regularly.

The liquid material should not be exposed to light at all, in order to avoid deterioration of the material's quality. Deviations from the explained manufacturing process can lead to altered mechanical properties of the FotoCast® material and color deviations.

Contains: 1,4 Butanediol dimethacrylate, Diphenyl(2,4,6-trimethylbenzoyl)phosphinoxide

Possible dangers:

Before use observe the Material Safety Data Sheet of FotoCast®.

Hazard statements

H317 May cause an allergic skin reaction.
H413 May cause long lasting harmful effects to aquatic life.

Precautionary statements

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P273 Avoid release to the environment.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.
P501 Dispose of contents/container in accordance with local/regional/national/international regulations.

Lot number / Expiry date: The lot number and the expiry date are indicated on each FotoCast® packaging. In case of claims please always indicate the lot number of the product. Do not use the product after the expiry date.

Dreve Otoplastik GmbH is not liable for any damages caused by improper application of the material.

2018-07-11
Dreve Otoplastik GmbH
Max-Planck-Str. 31 • 59423 Unna / Germany
Tel.: + 49 2303 8807-0 • Fax: + 49 2303 82909
E-Mail: otoplastik@dreve.de • Internet: www.dreve.com

Warning